

FD NeuroSilver™ Kit II

A rapid silver staining kit for the microscopic
detection of neuronal damage

User Manual
PK 301/301A, Version 2014-01

FOR IN VITRO RESEARCH USE ONLY
not for diagnostic or other uses



Quality & Excellence Since 1996

FD NeuroTechnologies
Consulting & Services, Inc.

© 1997-2014 All Rights Reserved

Contents

I.	Introduction	4
II.	Kit Contents	5
III.	Materials Required but Not Included	5
IV.	Safety and Handling Precautions	6
V.	Tissue Preparation	6
VI.	Staining Procedure	8
VII.	Evaluation of Results	11
VIII.	References using this kit	12

I. Introduction

FD NeuroSilver™ Kit II, an improved FD NeuroSilver™ Kit I, is designed for the detection of degenerating neurons in fixed tissue sections of the central nervous system from experimental animals. The principle of this kit is based on the findings that certain components of neurons undergoing degeneration, such as lysosomes, axons and terminals, become particularly argyrophilic. Under certain conditions, these cellular elements bind to silver ions with high affinity. Upon reduction, the silver ions form metallic grains that are visible under a light or an electron microscope.

FD NeuroSilver™ Kit has been widely used in animal studies under various experimental conditions (cf. page 11 for references using this kit). This kit has proven to be very specific and sensitive for the detection of degenerating neuronal somata, axons and terminals in both the brain and the spinal cord. It is particularly useful for the detection of small numbers of degenerating neurons that may not be demonstrable with routine histopathological techniques (for photo samples, please visit our website at www.fdneurotech.com).

FD NeuroSilver™ Kit has also been proven to be sensitive and reliable for the detection of amyloid plaques in the brains of transgenic mice. In addition, this kit may be used for demonstrating neurodegeneration and/or amyloid plaques in tissue sections that have been processed for immunohistochemistry (cf. *Nature Neuroscience* 3: 1301-1306, 2000).

II. Kit Contents

Store at room temperature

	<u>PK301A</u>	<u>PK301</u>
Solution A	250 ml	500 ml
Solution B	250 ml	500 ml
Solution C	250 ml x 2	500 ml x 2
Solution D	250 ml	500 ml
Solution E	1 ml	2 ml
Solution F	2 ml	3 ml
Solution G (10X)	250 ml	500 ml
Glass specimen retriever	2	2
Natural hair paintbrush	1	1
User Manual	1	1

III. Materials Required but Not Included

1. Double distilled or Milli-Q water
2. 6-well tissue culture plates
3. (Optional) Corning® Netwells® inserts (Corning Cat. #3480) and carrier kit (Corning Cat. #3521)
4. Histological supplies and equipment:
 - Microscope slides and coverslips
 - Staining jars
 - Xylenes
 - Permount®
 - A light microscope.

Corning® & Netwells® are registered trademarks of Corning, Inc.
Permount® is a registered trademark of Fisher Scientific.

IV. Safety and Handling Precautions

1. FD NeuroSilver™ Kit II is made for *in vitro* research use only and not for drug, diagnostic or other uses.
2. The kit contains reagents that may be harmful in contact with skin or by inhalation and fatal if ingested. Do not pipette by mouth. Avoid inhalation and contact with skin and eyes. In case of contact, wash immediately with generous amounts of water and seek medical advice. If swallowed, wash out mouth with water and immediately call a physician.
3. Perform experiment under a chemical hood. **Wear suitable protective clothing, gloves and eye/face protection while handling kit reagents.** Wash hands thoroughly after performing the experiment.

V. Tissue Preparation

Note

FD NeuroSilver™ Kit II has been proven to produce the best results in animal brain sections prepared according to the following procedure. However, it may be used in tissue specimens prepared differently (please contact FD NeuroTechnologies for more information).

1. Experimental animals should be deeply anesthetized and perfused via the ascending aorta with 0.1 M phosphate buffer (PB, pH 7.4) for 1 minute followed by 0.1 M PB containing 4% paraformalde-

hyde for 15 minutes (different species may require a different amount of perfusates, e.g. for a rat, 100 ml of saline and 500 ml of fixative should be sufficient). Upon removal, brains should be postfixed in the same fixative as that used for the perfusion overnight at 4°C. Subsequently, brains should be immersed in 0.1 M PB containing 20% or 30% sucrose for 72 hours at 4°C before freezing.

Note

If thick sections (e.g. above 60 µm) will be cut with a cryostat, we highly recommend the use of our tissue cryoprotection solution (Cat. #: PC102) instead of the sucrose solution.

2. Sections of 40 – 80 µm thickness may be cut either on a cryostat* or a vibratome (similar types of microtomes, such as a sliding microtome may be used alternatively). Sections should be collected in 0.1 M PB (pH 7.4) and be subsequently transferred into 0.1 M PB containing 4% paraformaldehyde. Store sections in 0.1 M PB containing 4% paraformaldehyde for at least 7 days at 4°C before processing with FD NeuroSilver™ kit II.

**To prevent tissue from ice crystal damage and to preserve the best possible cell morphology, tissue should be frozen rapidly before sectioning with a cryostat. For example, tissue may be rapidly frozen as described below: place tissue in a plastic spoon and slowly dip into iso-pentane pre-cooled with dry*

ice (for the best results, temperature should be kept below -70°C and the dipping should take about 1 min, the slower the better). After the tissue is completely immersed in iso-pentane, keep it in iso-pentane for a few seconds and then place it on dry ice for another minute to ensure that the tissue is well frozen. Do not let tissue thaw before sections are cut.

① Warning

Paraformaldehyde is harmful and toxic to humans if ingested or inhaled. The experiment with paraformaldehyde should be performed under a chemical hood with appropriate protection. Wear glasses and disposable gloves.

VI. Staining Procedure

The following instructions MUST be read before using this kit.

- *All containers (plastic preferred) to be used should be cleansed and rinsed with distilled water.*
- *Do not use metal implements whenever Solutions E is present.*
- *Prepare fresh working solutions and keep containers tightly closed at all times.*
- *Use 0.8 ml of solution for each 2 cm² of the section (e.g. approximately a rat brain coronal section through the dorsal hippocampus) to be stained at each step.*
- *Avoid bringing excess solution from one step to another. Very important!*

- *Do not let sections dry up anytime before mounting them on slides.*
- *Use of positive control sections is recommended for each experiment.*
- *The following procedure takes approximately 1 hour and should be performed at room temperature.*

1. Rinse sections in double distilled or Milli-Q water 2 times, 5 minutes each.
2. Transfer sections into a mixture containing equal volumes of **Solutions A** and **B**, 2 times, 10 minutes each.

 **Note**

- *The incubation plate must be covered during incubation to prevent splash and vaporization of the reagent.*
- *Sections should be immersed in the solution and incubated free floating with brisk shakes.*

3. Place sections in a mixture consisting of equal volumes of **Solutions A** and **B** with **Solution E** (one drop of Solution E for each 8 ml of a total volume of Solutions A and B) for 10 minutes.

e.g.	Solution A	4 ml
	Solution B	4 ml
	Solution E	1 drop

 **Note**

- *All 3 solutions must be well mixed and the resultant mixture should be totally clear without precipitates.*
- *The incubation plate must be covered during incubation to prevent splash and vaporization of the reagent.*
- *Sections should be immersed in the solution and incubated free floating with brisk shakes.*

4. Place sections in a mixture of **Solution C** and **Solution F** (one drop of Solution F for each 25 ml of Solution C), 2 times, 2 minutes each*.

** A total of 4 minutes of incubation time is satisfactory in most cases. However, variations in tissue processing may require that the duration of this step be lengthened (for decreasing the background stain) or shortened (for increasing the staining intensity). Note that prolonging the time of incubation will also decrease or possibly eliminate specific silver staining.*

■ Note

- *The incubation plate must be covered during incubation to prevent splash and vaporization of the reagent.*
- *Sections should be immersed in the solution and incubated free floating with brisk shakes.*

5. Place sections in a mixture of **Solution D** and **Solution F** (one drop of Solution F for each 25 ml of Solution D) for 5 minutes.

■ Note

- *The incubation plate must be covered during incubation to prevent splash and vaporization of the reagent.*
- *Sections should be immersed in the solution and incubated free floating with brisk shakes.*

6. Rinse sections in double distilled or Milli-Q water 2 times, 3 minutes each.
7. Place sections in 1x **Solution G** (10x concentrated solution provided, dilute with double distilled or Milli-Q water before use), 2 times, 5 minutes each.

Note

- *For the best results, keep sections in 1x Solution G at 4°C for 1-2 hours or longer if the background is high before mounting.*
- *Silver-stained sections should be protected from light whenever possible.*

8. Mount sections on pre-cleaned slides directly from 1x Solution G and let air-dry.

Note

For easy mounting, do not use adhesive microscope slides, including Superfrost Plus slides unless sections tend to fall off the slides during dehydration.

9. Clear sections in xylenes, 3 changes, 3 minutes each, and coverslip with Permount®.

Note

Do not dehydrate sections in ethanol as ethanol may cause loss of the staining.

VII. Evaluation of Results

Sections should be examined under a light microscope with either a brightfield or a darkfield condenser. Neurons undergoing degeneration are indicated by dense silver precipitates, appearing as black (brightfield) or golden (darkfield) grains, in their somata and/or processes. Photo samples are available at www.fdneurotech.com.

VIII. References using this kit:

1. Ahima RS, Bjorbaek C, Osei S, Flier JS (1999) Regulation of neuronal and glial proteins by leptin: implications for brain development. *Endocrinology* 140:2755-62.
2. Burns LH, Jin Z, Bowersox SS (1999) The neuroprotective effects of intrathecal administration of the selective N-type calcium channel blocker ziconotide in a rat model of spinal ischemia. *J. Vasc. Surg.* 30:334-43.
3. Betarbet R, Sherer TB, MacKenzie G, Garcia-Osuna M, Panov AV, Greenamyre JT (2000) Chronic systemic pesticide exposure reproduces features of Parkinson's disease. *Nature Neurosci.* 3:1301-6.
4. Kuhlmann AC, Guilarde TR (2000) Cellular and subcellular localization of peripheral benzodiazepine receptors after trimethyltin neurotoxicity. *J. Neurochem.* 74:1694-704.
5. Sternfeld M, Shoham S, Klein O, Flores-Flores C, Evron T, Idelson GH, Kitsberg D, Patrick JW, Soreq H (2000) Excess "read-through" acetylcholinesterase attenuates but the "synaptic" varian intensifies neurodegeneration correlates. *Proc. Natl. Acad. Sci. USA* 97:8647-52.
6. Ding Y, McAllister JP, Yao B, Yan N, Canady AI (2001) Axonal damage associated with enlargement of ventricles during hydrocephalus: a silver impregnation study. *Neurol. Res.* 23:581-7.
7. Ding Y, McAllister JP, Yao B, Yan N, Canady AI (2001) Neuron tolerance during hydrocephalus. *Neurosci.* 106:659-67.
8. Ding Y, Yao B, Lai Q, McAllister JP (2001) Impaired motor learning and diffuse axonal damage in motor and visual systems of the rat following traumatic brain injury. *Neurol. Res.* 23:193-202.
9. Freeman A, Ciliax B, Bakay R, Daley J, Miller RD, Keating G, Levey A, Rye D (2001) Nigrostriatal collaterals to thalamus degenerate in Parkinsonian animal models. *Ann. Neurol.* 50:321-9.
10. Liu ZP, Gastard M, Verina T, Bora S, Mouton PR, Koliatsos VE (2001) Estrogens modulate experimentally induced apoptosis of granule cells in the adult hippocampus. *J. Comp. Neurol.* 441:1-8.
11. Northington FJ, Ferriero DM, Martin LJ (2001) Neurodegeneration in the thalamus following neonatal hypoxia-ischemia is programmed cell death. *Dev. Neurosci.* 23:186-91.
12. Northington FJ, Ferriero DM, Graham EM, Traystman RJ, Martin LJ (2001) Early neurodegeneration after hypoxia-ischemia in neonatal rat is necrosis while delayed neuronal death is apoptosis. *Neurobiol. Dis.* 8:207-19.
13. Shoham S, Javitt DC, Heresco-Levy U (2001) Chronic high-dose glycine nutrition: effects on rat brain cell morphology. *Biol. Psychiat.* 49:876-85.
14. Zito MA, Koenneke LA, McAuliffe MJ, McNally B, van Rooijen N, Heyes MP (2001) Depletion of systemic macrophages by liposome-encapsulated clodronate attenuates striatal macrophage invasion and neurodegeneration following local endotoxin infusion in gerbils. *Brain Res.* 892:13-26.
15. Abdel-Rahman A, Shetty AK, Abou-Donia MB (2002) Disruption of the blood-brain barrier and neuronal cell death in cingulate cortex, dentate gyrus, thalamus, and hypothalamus in a rat model of gulf-war syndrome. *Neurobiol. Dis.* 10:306-26.
16. Bianco CL, Ridet JL, Schneide BL, Déglon N, Aebscher P (2002) α -synucleinopathy and selective dopaminergic neuron loss in a rat lentiviral-based model of Parkinson's disease. *Proc. Natl. Acad. Sci. USA* 99:10813-8.
17. Ding Y, Yao B, Zhou Y, Park H, McAllister JP, Diaz FG (2002) Prereperfusion flushing of ischemic territory: a therapeutic study in which histological and behavioral assessments were used to measure ischemia-reperfusion injury in rats with stroke. *J. Neurosurg.* 96:310-9.
18. Gordon MN, Holcomb LA, Jantzen PT, DiCarlo G, Wilcock D, Boyett KW, Connor K, Melachrino J, O'Callaghan JP, Morgan D (2002) Time course of the development of Alzheimer-like pathology in the doubly transgenic PS1+APP mouse. *Exp. Neurol.* 173:183-95.
19. Guo-Ross SX, Clark S, Montoya DAC, Jones KH, Obernier J, Shetty AK, White AM, Bluszta JK, Wilson WA, Swartzwelder HS (2002) Prenatal choline supplementation protects against postnatal neurotoxicity. *J. Neurosci.* 22:RC195:1-6.
20. McCormack AL, Thiruchelvam M, Manning-Bog AB, Thiffault C, Langston JW, Cory-Slechta DA, Di Monte DA (2002) Environmental risk factors and Parkinson's disease: selective degeneration of nigral dopaminergic neurons caused by the herbicide paraquat. *Neurobiol. Dis.* 10:119-27.

21. Montoliu C, Humet M, Canales J-J, Burda J, Planells-Cases R, Sánchez-Baeza F, Carbonell T, Pérez-Payá E, Messeguer A, Ferrer-Montiel A, Felipo V (2002) Prevention of in vivo excitotoxicity by a family of trialkylglycines, a novel class of neuroprotectants. *J. Pharmacol. Exp. Ther.* 301:29-36.
22. Tong W, Igarashi T, Ferriero DM, Noble LJ (2002) Traumatic brain injury in the immature mouse brain: characterization of regional vulnerability. *Exp. Neurol.* 176:105-16.
23. Lauwers E, Debyser Z, Van Dorpe J, De Strooper B, Nuttin B, Baekelandt V (2003) Neuropathology and neurodegeneration in rodent brain induced by lentiviral vector-mediated overexpression of alpha-synuclein. *Brain Pathol.* 13:364-72.
24. Manning-Bog AB, McCormack AL, Purisai MG, Bolin LM, Di Monte A (2003) α -synuclein overexpression protects against paraquat-induced neurodegeneration. *J. Neurosci.* 23:3095-9.
25. Sherer TB, Kim JH, Betarbet R, Greenamyre JT (2003) Subcutaneous rotenone exposure causes highly selective dopaminergic degeneration and α -synuclein aggregation. *Exp. Neurol.* 179:9-16.
26. Shoham S, Bejar C, Kovalev E, Weinstock M (2003) Intracerebroventricular injection of streptozotocin causes neurotoxicity to myelin that contributes to spatial memory deficits in rats. *Exp. Neurol.* 184: 1043-52.
27. Abdel-Rahman A, Abou-Donia SM, El-Masry EM, Shetty AK, Abou-Donia MB (2004) Stress and combined exposure to low doses of pyridostigmine bromide, deet, and permethrin produce neurochemical and neuropathological alterations in cerebral cortex, hippocampus, and cerebellum. *J. Toxicol. Env. Heal. A.* 67:163-92.
28. Lo Bianco C, Schneider BL, Bauer M, Sajadi A, Brice A, Iwatsubo T, Aebsicher P (2004) Lentiviral vector delivery of parkin prevents dopaminergic degeneration in an α -synuclein rat model of Parkinson's disease. *Proc. Natl. Acad. Sci. USA* 101:17510-5.
29. Paschen W, Yatsiv I, Shoham S, Shohami E (2004) Brain trauma induces X-box protein 1 processing indicative of activation of the endoplasmic reticulum unfolded protein response. *J. Neurochem.* 88:983-92.
30. Shoham S, Youdim MBH (2004) Nutritional iron deprivation attenuates kainate-induced neurotoxicity in rats: implications for involvement of iron in neurodegeneration. *Ann. N.Y. Acad. Sci.* 1012:24-114.
31. Liang ZQ, Wang XX, Wang Y, Chuang DM, DiFiglia M, Chase TN, Qin ZH (2005) Susceptibility of striatal neurons to excitotoxic injury correlates with basal levels of Bcl-2 and the induction of P53 and c-Myc immunoreactivity. *Neurobiol. Dis.* 20:562-73.
32. Williams AJ, Hartings JA, Lu XCM, Rolli ML, Dave JR, Tortella FC (2005) Characterization of a new rat model of penetrating ballistic brain injury. *J. Neurotrauma* 22:313-31.
33. Adhami F, Liao GH, Morozov YM, Schloemer A, Schmitzhorst VJ, Lorenz JN, Dunn RS, Vorhees CV, Wills-Karp M, Degen JL, Davis RJ, Mizushima N, Rakic P, Dardzinski BJ, Holland SK, Sharp FR, Kuan CY (2006) Cerebral ischemia-hypoxia induces intravascular coagulation and autophagy. *Am. J. Pathol.* 169:566-83.
34. Cui L, Jeong H, Borovecki F, Parkhurst CN, Tanese N, Krainc D (2006) Transcriptional repression of PGC-1 α by mutant huntingtin leads to mitochondrial dysfunction and neurodegeneration. *Cell* 127:59-69.
35. Dow G, Bauman R, Caridha D, Cabezas M, Du F, Gomez-Lobo R, Park M, Smith K, Cannard K (2006) Mefloquine induces dose-related neurological effects in a rat model. *Antimicrob. Agents Ch.* 50:1045-53.
36. Luo J, Lin AH, Masliah E, Wyss-Coray T (2006) Bioluminescence imaging of Smad signaling in living mice shows correlation with excitotoxic neurodegeneration. *Proc. Natl. Acad. Sci. USA* 103:18326-31.
37. Martin LJ, Pan Y, Price AC, Sterling W, Copeland NG, Jenkins NA, Price DL, Lee MK (2006) Parkinson's disease α -synuclein transgenic mice develop neuronal mitochondrial degeneration and cell death. *J. Neurosci.* 26:41-50.
38. Rao MS, Hattiangady B, Reddy DS, Shetty AK (2006) Hippocampal neurodegeneration, spontaneous seizures, and mossy fiber sprouting in the F344 rat model of temporal lobe epilepsy. *J. Neurosci. Res.* 83:1088-105.
39. Saydoff JA, Garcia RAG, Browne SE, Liu LS, Sheng J, Brenneman D, Hu Z, Cardin S, Gonzalez A, von Borstel RW, Gregorio J, Burr H, Flint Beal M (2006) Oral uridine pro-drug PN401 is neuroprotective in the R6/2 and N171-82Q mouse models of Huntington's disease. *Neurobiol. Dis.* 24:455-65.
40. Theodore S, Stolberg S, Cass WA, Maragos WF (2006) Human immunodeficiency virus-1 protein Tat and methamphetamine interactions. *Ann. N.Y. Acad. Sci.* 1074:178-90.
41. Williams AJ, Hartings JA, Lu XCM, Rolli ML, Tortella FC (2006) Penetrating ballistic-like brain injury in the rat: differential time courses of hemorrhage, cell death, inflammation, and remote degeneration. *J. Neurotrauma* 23:1828-46.

42. Yan X, Benson JM, Gomez AP, Baden DG, Murray TF (2006) Brevetoxin-induced neural insult in the retrosplenial cortex of mouse brain. *Inhal. Toxicol.* 18:1109-16.
43. Yazdani U, German DC, Liang C-L, Manzino L, Sonsalla PK, Zeevalk GD (2006) Rat model of Parkinson's disease: chronic central delivery of 1-methyl-4-phenylpyridinium (MPP⁺). *Exp. Neurol.* 200:172-83.
44. Zhang Y, Leavitt BR, van Raamsdonk JM, Dragatsis I, Goldowitz D, MacDonald ME, Hayden MR, Friedlander RM (2006) Huntingtin inhibits caspase-3 activation. *EMBO J.* 25:5896-906.
45. Bliss RM, Hoffman T, Finchbone V, Zeitouni T, Mansouri B, Strahlendorf J (2007) AMPA-induced excitotoxicity stimulates calpain activity and is exacerbated by tumor necrosis factor- α (TNF- α) in the CA1 and CA3 regions of the hippocampus in an acute slice model. *FASEB J.* 21:9432.
46. Caudle WM, Richardson JR, Wang MZ, Taylor TN, Guillot TS, McCormack AL, Colebrooke RE, Di Monte DA, Emson PC, Miller GW (2007) Reduced vesicular storage of dopamine causes progressive nigrostriatal neurodegeneration. *J. Neurosci.* 27:8138-48.
47. Kauk D, Rajagopalan S, Chinta S, Kumar J, Monto DD, Cherny RA, Andersen JK (2007) Chronic ferritin expression within murine dopaminergic midbrain neurons results in a progressive age-related neurodegeneration. *Brain Res.* 1140:188-94.
48. Liang ZQ, Li YL, Zhao XL, Han R, Wang XX, Wang Y, Chase TN, Bennett MC, Qin ZH (2007) NF-kappaB contributes to 6-hydroxydopamine-induced apoptosis of nigral dopaminergic neurons through p53. *Brain Res.* 1145:190-203.
49. Parthasarathy KS, Cheng YC, McAllister JP 2nd, Shen Y, Li J, Deren K, Haacke EM, Auner GW (2007) Biocompatibilities of sapphire and borosilicate glass as cortical neuroprostheses. *Magn. Reson. Imaging* 25(9):1333-40.
50. Setsuie R, Wang YL, Mochizuki H, Osaka H, Hayakawa H, Ichihara N, Li H, Furuta A, Sano Y, Sun YJ, Kwon J, Kabuta T, Yoshimini K, Aoki S, Mizuno Y, Noda M, Wada K (2007) Dopaminergic neuronal loss in transgenic mice expressing the Parkinson's disease-associated UCH-L1 193M mutant. *Neurochem. Int.* 50:119-29.
51. Wang A, Liang X, McAllister II JP, Li J, Brabant K, Black C, Finlayson P, Cao T, Tang H, Salley SO, Auner GW, Simon Ng KY (2007) Stability of and inflammatory response to silicon coated with a fluoroalkyl self-assembled monolayer in the central nervous system. *J. Biomed. Mater. Res. A* 81:363-72.
52. Williams AJ, Wei HH, Dave JR, Tortella FC (2007) Acute and delayed neuroinflammatory response following experimental penetrating ballistic brain injury in the rat. *J. Neuroinflammation* 4:17.
53. Xu W, Wong TP, Chery N, Gaertner T, Wang YT, Baudry M (2007) Calpain-mediated mGluR1 α truncation: a key step in excitotoxicity. *Neuron* 53:399-412.
54. Berman RF, Pessah IN, Mouton PR, Mav D, Harry J (2008) Low-level neonatal thimerosal exposure: further evaluation of altered neurotoxic potential in SJL mice. *Toxicol. Sci.* 101 (2):294-309.
55. Calingasan NY, Ho DJ, Wille EJ, Campagna MV, Ruan J, Dumont M, Yang L, Shi Q, Gibson GE, Beal MF (2008) Influence of mitochondrial enzyme deficiency on adult neurogenesis in mouse models of neurodegenerative diseases. *Neurosci.* 153(4):986-96.
56. Chen L, Ding Y, Cagniard B, Van Laar AD, Mortimer A, Chi W, Hastings TG, Kang UJ, Zhuang X (2008) Unregulated cytosolic dopamine causes neurodegeneration associated with oxidative stress in mice. *J. Neurosci.* 28(2):425-33.
57. Chen ZL, Yu H, Yu WM, Pawlak R, Strickland S (2008) Proteolytic fragments of laminin promote excitotoxic neurodegeneration by up-regulation of the KA1 subunit of the kainate receptor. *J. Cell Biol.* 183(7):1299-313.
58. Dombrowski SM, Deshpande A, Dingwall C, Leichliter A, Leibson Z, Luciano MG (2008) Chronic hydrocephalus-induced hypoxia: increased expression of VEGFR-2+ and blood vessel density in hippocampus. *Neurosci.* 152(2):346-59.
59. Eid T, Ghosh A, Wang Y, Beckström H, Zaveri HP, Lee TS, Lai JC, Mathankar-Phatak GH, de Lanerolle NC (2008) Recurrent seizures and brain pathology after inhibition of glutamine synthetase in the hippocampus in rats. *Brain* 131(Pt 8):2061-70.
60. Flanagan-Cato LM, Fluharty SJ, Weinreb EB, LaBelle DR (2008) Food restriction alters neuronal morphology in the hypothalamic ventromedial nucleus of male rats. *Endocrinology* 149(1):93-9.
61. GuilarTE TR, Burton NC, Verina T, Prabhu VV, Becker KG, Syversen T, Schneider JS (2008) Increased APLP1 expression and neurodegeneration in the frontal cortex of manganese-exposed non-human primates. *J. Neurochem.* 105(5):1948-59.
62. Guillot TS, Shepherd KR, Richardson JR, Wang MZ, Li Y, Emson PC, Miller GW (2008) Reduced vesicular storage of dopamine exacerbates methamphetamine-induced neurodegeneration and astrogliosis. *J. Neurochem.* 106(5):2205-17.

63. Jackson N, Muthuswamy J (2008) Artificial dural sealant that allows multiple penetrations of implantable brain probes. *J. Neurosci. Methods.* 171(1):147-52.
64. Lehrmann E, Guidetti P, Löve A, Williamson J, Bertram EH, Schwarcz R (2008) Glial activation precedes seizures and hippocampal neurodegeneration in measles virus-infected mice. *Epilepsia* 49 Suppl 2:13-23.
65. Liu L, Orozco JJ, Planel E, Wen Y, Bretteville A, Krishnamurthy P, Wang L, Herman M, Figueira H, Yu WH, Arancio O, Duff K (2008) A transgenic rat that develops Alzheimer's disease-like amyloid pathology, deficits in synaptic plasticity and cognitive impairment. *Neurobiol. Dis.* 31(1):46-57.
66. Mallajosyula JK, Kau D, Chinta SJ, Rajagopalan S, Rane A, Nicholls DG, Di Monte DA, Macarthur H, Andersen JK (2008) MAO-B elevation in mouse brain astrocytes results in Parkinson's pathology. *PLoS One* 3(2):e1616.
67. Masuda N, Peng Q, Li Q, Jiang M, Liang Y, Wang X, Zhao M, Wang W, Ross CA, Duan W (2008) Tiagabine is neuroprotective in the N171-82Q and R6/2 mouse models of Huntington's disease. *Neurobiol. Dis.* 30(3):293-302.
68. Moussa CE, Rusnak M, Hailu A, Sidhu A, Fricke ST (2008) Alterations of striatal glutamate transmission in rotenone-treated mice: MRI/MRS *in vivo* studies. *Exp. Neurol.* 209 (1):224-33.
69. Wetzel MK, Naska S, Laliberté CL, Rymer VV, Fujitani M, Biernaskie JA, Cole CJ, Lerch JP, Spring S, Wang SH, Frankland PW, Henkelman RM, Josselyn SA, Sadikot AF, Miller FD, Kaplan DR (2008) p73 regulates neurodegeneration and phospho-tau accumulation during aging and Alzheimer's disease. *Neuro* 59(5):708-21.
70. Bauman RA, Ling G, Tong L, Januszkiewicz A, Agoston D, de Lanerolle N, Kim Y, Ritzel D, Bell R, Ecklund J, Armonda R, Bandak P, Parks S (2009) An introductory characterization of a combat-casualty-care relevant swine model of closed head injury resulting from exposure to explosive blast. *J. Neurotrauma* 26(6):841-60.
71. Burns MP, Zhang L, Rebeck GW, Querfurth HW, Moussa CE (2009) Parkin promotes intracellular Abeta1-42 clearance. *Hum. Mol. Genet.* 18(17):3206-16.
72. Choi DY, Liu M, Hunter RL, Cass WA, Pandya JD, Sullivan PG, Shin EJ, Kim HC, Gash DM, Bing G (2009) Striatal neuroinflammation promotes Parkinsonism in rats. *PLoS One* 4(5):e5482.
73. Deshpande A, Dombrowski SM, Leichliter A, Krajeir N, Zingales N, Inoue M, Schenck S, Fukamachi K, Luciana MG (2009) Dissociation between vascular endothelial growth factor receptor-2 and blood vessel density in the caudate nucleus after chronic hydrocephalus. *J. Cerebr. Blood F. Met.* 29(11):1806-15.
74. Han SR, Shin C, Park S, Rhyu S, Park J, Kim YI (2009) Differential expression of activating transcription factor-2 and c-Jun in the immature and adult rat hippocampus following lithium-pilocarpine induced status epilepticus. *Yonsei. Med. J.* 50(2):200-5.
75. Jinde S, Belforte JE, Yamamoto J, Wilson MA, Tonegawa S, Nakazawa K (2009) Lack of kainic acid-induced gamma oscillations predicts subsequent CA1 excitotoxic cell death. *Eur. J. Neurosci.* 30(6):1036-55.
76. Long JB, Bentley TL, Wessner KA, Cerone C, Sweeney S, Bauman RA (2009) Blast overpressure in rats: recreating a battlefield injury in the laboratory. *J. Neurotrauma* 26 (6):827-40.
77. Pernot F, Carpentier P, Baille V, Testylier G, Beaufort C, Foquin A, Filliat P, Liscia P, Coutan M, Piérard C, Béracochea D, Dorandeu F (2009) Intrahippocampal cholinesterase inhibition induces epileptogenesis in mice without evidence of neurodegenerative events. *Neurosci.* 162(4):1351-65.
78. Rajasekaran K, Sun C, Bertram EH (2009) Altered pharmacology and GABA-A receptor subunit expression in dorsal midline thalamic neurons in limbic epilepsy. *Neurobiol. Dis.* 33(1):119-32.
79. Shear DA, Williams AJ, Sharow K, Lu XC, Tortella FC (2009) Neuroprotective profile of dextromethorphan in an experimental model of penetrating ballistic-like brain injury. *Pharmacol. Biochem. Behav.* 94(1):56-62.
80. Sloan DM, Bertram EH (2009) Changes in midline thalamic recruiting responses in the prefrontal cortex of the rat during the development of chronic limbic seizures. *Epilepsia* 50 (3):556-65.
81. Vig PJ, Shao Q, Subramony SH, Lopez ME, Safaya E (2009) Bergmann glial S100B activates myo-inositol monophosphatase 1 and co-localizes to purkinje cell vacuoles in SCA1 transgenic mice. *Cerebellum* 8(3):231-44.
82. Valor LM, Jancic D, Lujan R, Barco A (2010) Ultrastructural and transcriptional profiling of neuropathological misregulation of CREB function. *Cell Death Differ.* 17(10):1636-44.
83. De Araujo Furtado M, Lumley LA, Robison C, Tong LC, Lichtenstein S, Yourick DL (2010) Spontaneous recurrent seizures after status epilepticus induced by soman in Sprague

- Dawley rats. *Epilepsia* 51(8):1503-10.
84. Khandelwal PJ, Dumanis SB, Feng LR, Maguire-Zeiss K, Rebeck GW, Lashuel HA, Moussa CEH (2010) Parkinson-related parkin reduces alpha-Synuclein phosphorylation in a gene transfer model. *Mol. Neurodegeneration* 5(47):1-13.
 85. Kocherhans S, Madhusudan A, Doehner J, Breu KS, Nitsch RM, Fritschy JM, Knuesel I (2010) Reduced Reelin expression accelerates amyloid-beta plaque formation and tau pathology in transgenic Alzheimer's disease mice. *J. Neurosci.* 30(27):9228-40.
 86. Kunimoto S, Nakamura S, Wada K, Inoue T (2010) Chronic stress-mutated presenilin 1 gene interaction perturbs neurogenesis and accelerates neurodegeneration. *Exp. Neurol.* 221(1):175-85.
 87. Liu M, Choi D-Y, Hunter RL, Pandya JD, Cass WA, Sullivan PG, Kim H-C, Gash DM, Bing G (2010) Trichloroethylene induces dopaminergic neurodegeneration in Fisher 344 rats. *J. Neurochem.* 112(3):773-83.
 88. Obeid M, Frank J, Medina M, Finckbone V, Bliss R, Bista B, Majmudar S, Hurst D, Strahlendorf H, Strahlendorf J (2010) Neuroprotective effects of leptin following kainic acid-induced status epilepticus. *Epilepsy Behav.* 19(3):278-83.
 89. Rebeck GW, Hoe H-S, Moussa CEH (2010) β -Amyloid Gene Transfer Model Exhibits Intraneuronal Amyloid, Gliosis, Tau Phosphorylation, and Neuronal Loss. *J. Biol. Chem.* 285:7440-6.
 90. Siddiqui A, Mallajosyula JK, Rane A, Andersen JK (2010) Ability to delay neuropathological events associated with astrocytic MAO- B increase in a Parkinsonian mouse model: implications for early intervention on disease progression. *Neurobiol. Dis.* 40(2):444-8.
 91. Smith WW, Liu Z, Liang Y, Masuda N, Swing DA, Jenkins NA, Copeland NG, Troncoso JC, Pletnikov M, Dawson TM, Martin LJ, Moran TH, Lee MK, Borchelt DR, Ross CA (2010) Synphilin-1 attenuates neuronal degeneration in the A53T alpha-synuclein transgenic mouse model. *Hum. Mol. Genet.* 19(11):2087-98.
 92. Wong M, Martin LJ (2010) Skeletal muscle-restricted expression of human SOD1 causes motor neuron degeneration in transgenic mice. *Hum. Mol. Genet.* 19(11):2284-302.
 93. Xu YF, Gendron TF, Zhang YJ, Lin WL, D'Alton S, Sheng H, Casey MC, Tong J, Knight J, Yu X, Rademakers R, Boylan K, Hutton M, McGowan E, Dickson DW, Lewis J, Petrucci L (2010) Wild-type human TDP-43 expression causes TDP-43 phosphorylation, mitochondrial aggregation, motor deficits, and early mortality in transgenic mice. *J. Neurosci.* 30(32):10851-9.
 94. Zhang J, Peng Q, Li Q, Jahanshad N, Hou Z, Jiang M, Masuda N, Langbehn DR, Miller MI, Mori S, Ross CA, Duan W (2010) Longitudinal characterization of brain atrophy of a Huntington's disease mouse model by automated morphological analyses of magnetic resonance images. *Neuroimage* 49(3):2340-51.
 95. Zhou H, Huang C, Chen H, Wang D, Landel CP, Xia PY, Bowser R, Liu YJ, Xia XG (2010) Transgenic rat model of neurodegeneration caused by mutation in the TDP gene. *PLoS Genet.* 6(3):e1000887.
 96. Bliss RM, Finckbone VL, Trice J, Strahlendorf H, Strahlendorf J (2011) Tumor necrosis factor- α (TNF- α) augments AMPA-induced Purkinje neuron toxicity. *Brain Res.* 1386:1-14.
 97. Davidsson J, Risling M (2011) A new model to produce sagittal plane rotational induced diffuse axonal injuries. *Front. Neurol.* 2:41.
 98. Khandelwal PJ, Herman AM, Hoe HS, Rebeck GW, Moussa CE (2011) Parkin mediates beclin-dependent autophagic clearance of defective mitochondria and ubiquitinated Abeta in AD models. *Hum. Mol. Genet.* 20(11):2091-102.
 99. Koliatios VE, Cernak I, Xu L, Song Y, Savonenko A, Crain BJ, Eberhart CG, Frangakis CE, Melnikova T, Kim H, Lee D (2011) A mouse model of blast injury to brain: initial pathological, neuropathological, and behavioral characterization. *J. Neuropathol. Exp. Neurol.* 70(5):399-416.
 100. Langston JL, Wright LK, Connis N, Lumley LA (2011) Characterizing the behavioral effects of nerve agent-induced seizure activity in rats: Increased startle reactivity and perseverative behavior. *Pharmacol Biochem Behav.* 100(3), 382-391.
 101. Moffett MC, Schultz MK, Schwartz JE, Stone MF, Lumley LA (2011) Impaired auditory and contextual fear conditioning in soman-exposed rats. *Pharmacol Biochem Behav.* 98 (1):120-9.
 102. Shitaka Y, Tran HT, Bennett RE, Sanchez L, Levy MA, Dikranian K, Brody DL (2011) Repetitive closed-skull traumatic brain injury in mice causes persistent multifocal axonal injury and microglial reactivity. *J Neuropathol Exp. Neurol.* 70(7):551-67.
 103. Wang Y, Pan Y, Price A, Martin LJ (2011) Generation and characterization of transgenic mice expressing mitochondrial targeted red fluorescent protein selectively in neurons: modeling mitochondriopathy in excitotoxicity and amyotrophic lateral sclerosis. *Molecular*

Neurodegeneration 6:75.

104. Wang Y, Wei Y, Oguntayo S, Wilkins W, Arun P, Valiyaveettil M, Song J, Long JB, Nambiar MP (2011) Tightly coupled repetitive blast-induced traumatic brain injury: development and characterization in mice. *J. Neurotrauma* 28(10):2171-83.
105. Acosta MT, Munashings J, Zhang L, Guerron DA, Vortmeyer A, Theodore WH (2012) Isolated seizures in rats do not cause neuronal injury. *Acta Neurol. Scand.* 125(1):30-37.
106. Bennett RE, Mac Donald CL, Brody DL (2012) Diffusion tensor imaging detects axonal injury in a mouse model of repetitive closed-skull traumatic brain injury. *Neurosci. Lett.* 513(2):160-5.
107. Brunjes PC (2012) The mouse olfactory peduncle. 2. The anterior limb of the anterior commissure. *Front. Neuroanat.* 6:51.
108. Doehner J, Genoud C, Imhof C, Krstic D, Knuessl I (2012) Extrusion of misfolded and aggregated proteins—a protective strategy of aging neurons? *Eur. J. Neurosci.* 35(12):1938-50.
109. Fünfschilling U, Jockusch WJ, Sivakumar N, Möbius W, Corthals K, Li S, Quintes S, Kim Y, Schaap IA, Rhee JS, Nave KA, Saher G (2012) Critical time window of neuronal cholesterol synthesis during neurite outgrowth. *J. Neurosci.* 32(22):7632-45.
110. Jiang Y, Brody DL (2012) Administration of COG1410 Reduces Axonal Amyloid Precursor Protein Immunoreactivity and Microglial Activation after Controlled Cortical Impact in Mice. *J. Neurotrauma* 29(13):2332-41.
111. Krstic D, Madhusudan A, Doehner J, Vogel P, Notter T, Imhof C, Manalastas A, Hilfiker M, Pfister S, Schwerdel C, Riether C, Meyer U, Knuessl I (2012) Systemic immune challenges trigger and drive Alzheimer-like neuropathology in mice. *J. Neuroinflammation*. 9 (1):151.
112. Langston JL, Wright LK, Connis N, Lumley LA (2012) Characterizing the behavioral effects of nerve agent-induced seizure activity in rats: Increased startle reactivity and perseverative behavior. *Pharmacol. Biochem. Behav.* 100(3):382-391.
113. Plantman S, Ng KC, Lu J, Davidsson J, Risling M (2012) Characterization of a Novel Rat Model of Penetrating Traumatic Brain Injury. *J. Neurotrauma* 29(6):1219-32.
114. Schultz MK, Wright LKM, Stone MF, Schwartz JE, Kelley NR, Moffett MC, Lee RB, Lumley LA (2012) The anticholinergic and antiglutamatergic drug caramiphen reduces seizure duration in soman-exposed rats: Synergies with the benzodiazepine diazepam. *Toxicol. App. Pharmacol.* 259(3):376-86.
115. Svetlov SI, Prima V, Glushakova O, Svetlov A, Kirk DR, Gutierrez H, Serebryany VL, Curley KC, Wang KK, Hayes RL (2012) Neuro-glial and systemic mechanisms of pathological responses in rat models of primary blast overpressure compared to "composite" blast. *Front. Neurol.* 3:15. PMID: 22403567.
116. Vignaux G, Chabbert C, Gaboyard-Niay S, Travo C, Machado ML, Denise P, Comoz F, Hitier M, Landemore G, Philoxene B, Besnard S (2012) Evaluation of the chemical model of vestibular lesions induced by arsanilate in rats. *Toxicol. App. Pharmacol.* 258(1):61-71.
117. Yoon SO, Park DJ, Ryu JC, Ozer HG, Tep C, Shin YJ, Lim TH, Pastorino L, Kunwar AJ, Walton JC, Nagahara AH, Lu KP, Nelson RJ, Tuszyński MH, Huang K (2012) JNK3 perpetuates metabolic stress induced by Aβ peptides. *Neuron* 75(5):824-37.
118. Feeney EJ, Stephenson D, Kleiman R, Bove S, Cron C, Moody L, Robinson M, Ramirez JJ (2013) Immunohistochemical characterization of axonal sprouting in mice. *Restorative Neurology and Neuroscience* 31(5):517-531.
119. Foxworthly WA, Clemo HR, Meredith MA (2013) Laminar and connectional organization of a multisensory cortex. *J. Comp. Neurol.* 521(8):1867-90.
120. Hylin MJ, Orsi SA, Zhao J, Bockhorst K H, Perez AI, Moore AN, Dash PK (2013) Behavioral and histopathological alterations resulting from mild fluid percussion injury. *J. Neurotrauma* 30(9):702-715
121. Iglesias-Bregna D, Hanak S, Ji Z, Petty M, Liu L, Zhang D, McMonagle-Strucko K (2013) Effects of Prophylactic and Therapeutic Teriflunomide in Transcranial Magnetic Stimulation-Induced Motor-Evoked Potentials in the Dark Agouti Rat Model of Experimental Autoimmune Encephalomyelitis. *J. Pharmacol. Exp. Ther.* 347(1):203-11.
122. Kraft AW, Hu X, Yoon H, Yan P, Xiao Q, Wang Y, Gil SC, Brown J, Wilhelmsson U, Restivo JL, Cirrito JR, Holtzman DM, Kim J, Pekny M, Lee JM (2013) Attenuating astrocyte activation accelerates plaque pathogenesis in APP/PS1 mice. *FASEB J.* 27(1):187-98.
123. Namjoshi DR, Martin G, Donkin J, Wilkinson A, Stukas S, Fan J, Carr M, Tabarestani S, Wuerth K, Hancock RE, Wellington CL (2013) The Liver X Receptor Agonist GW3965 Improves Recovery from Mild Repetitive Traumatic Brain Injury in Mice Partly through Apolipoprotein E. *PLoS One* 8(1): e53529.
124. Tong J, Huang C, Bi F, Wu Q, Huang B, Liu XH, Li F, Zhou HX, Xia XG (2013) Express-

- sion of ALS-linked TDP-43 mutant in astrocytes causes non-cell-autonomous motor neuron death in rats. **EMBO J.** 32(13):1917-26.
125. Baizer JS (2014) Unique features of the human brainstem and cerebellum. **Frontiers in Human Neuroscience** 8:202.
 126. Bennett RE, Brody DL (2014) Acute Reduction of Microglia Does Not Alter Axonal Injury In a Mouse Model of Repetitive Concussive Traumatic Brain Injury. **J. Neurotrauma** 31(19):1647-63.
 127. Byrns CN, Pitts MW, Gilman CA, Hashimoto AC, Berry MJ (2014) Mice Lacking Seleno-protein P and Selenocysteine Lyase Exhibit Severe Neurological Dysfunction, Neurodegeneration, and Audiogenic Seizures. **J. Biological Chemistry** 289(14):9662-74.
 128. Cachón-González MB, Wang SZ, Ziegler R, Cheng SH, Cox TM (2014) Reversibility of neuropathology in Tay-Sachs-related diseases. **Human Molecular Genetics** 23(3):730-48.
 129. Chen Z, Shin D, Chen S, Mikhail K, Hadass O, Tomlison BN, Korkin D, Shyu CR, Cui JK, Anthony DC, Gu Z (2014) Histological Quantitation of Brain Injury Using Whole Slide Imaging: A Pilot Validation Study in Mice. **PLoS One** 9(3):e92133.
 130. Hebron M, Chen W, Miessau MJ, Lonskaya I, Moussa CEH (2014) Parkin reverses TDP-43-induced cell death and failure of amino acid homeostasis. **J. Neurochem.** 129(2):350-61.
 131. Heldt SA, Elberger AJ, Deng Y, Guley NH, Del Mar N, Rogers J, Choi GW, Ferrell J, Rex TS, Honig MG, Reiner A (2014) A novel closed-head model of mild traumatic brain injury caused by primary overpressure blast to the cranium produces sustained emotional deficits in mice. **Frontiers in Neurology** 5:47.
 132. Kovacs SK, Leonessa F, Ling GS (2014) Blast TBI models, neuropathology, and implications for seizure risk. **Frontiers in Neurology** 5:47.
 133. Larivière R, Gaudet R, Gentil BJ, Girard M, Conte TC, Minotti S, Leclerc-Desauniers I K, Gehring K, McKinney RA, Shoubridge EA, McPherson PS, Durham HD, Brais B (2014) Sacs knockout mice present pathophysiological defects underlying autosomal recessive spastic ataxia of Charlevoix-Saguenay. **Human Molecular Genetics** ddu491.
 134. Martin LJ, Semenkova S, Hanford A, Wong M (2014) The mitochondrial permeability transition pore regulates Parkinson's disease development in mutant α -synuclein transgenic mice. **Neurobiology of Aging** 35(5):1132-52.
 135. Potter MC, Wozniak KM, Calizzot N, Slusher BS (2014) Glutamate Carboxypeptidase II Inhibition Behaviorally and Physiologically Improves Pyridoxine-Induced Neuropathy in Rats. **PLoS One** 9(9):e102936.
 136. Rozas de O'Laughlin N, Redell JB, Hill J, McKenna J, Moore AN, Gambello MJ, Dash PK (2014) Genetic activation of m TORC1 signaling worsens neurocognitive outcome after traumatic brain injury. **J. Neurotrauma** PMC3988378.
 137. Schultz MK, Wright LKM, de Araujo Furtado M, Stone MF, Moffett MC, Kelley NR, Bourne AR, Lumeh WZ, Schultz CR, Schwartz JE, Lumley LA (2014) Caramiphen edisylate as adjunct to standard therapy attenuates soman-induced seizures and cognitive deficits in rats. **Neurotoxicology and Teratology** 44:89-104.
 138. Tsika E, Glauer L, Moser R, Fiser A, Daniel G, Sheerin UM, Lees A, Troncoso JC, Lewis PA, Bandopadhyay R, Schneider BL, Moore DJ (2014) Parkinson's disease-linked mutations in VPS35 induce dopaminergic neurodegeneration. **Human Molecular Genetics** ddu178.
 139. Tsika E, Kannan M, Foo CSY, Dikeman D, Glauer L, Gellhaar S, Galter D, Knott GW, Dawson TM, Dawson VL, Moore DJ (2014) Conditional Expression of Parkinson's Disease-Related R1441C LRRK2 in Midbrain Dopaminergic Neurons of Mice Causes Nuclear Abnormalities without Neurodegeneration. **Neurobiol. Dis.** 71:345-58.
 140. Valiyaveettil M, Alamneh YA, Wang Y, Arun P, Oguntayo S, Wei YL, Long JB, Nambiar MP (2014) Cytoskeletal protein α -II spectrin degradation in the brain of repeated blast exposed mice. **Brain Res.** 1549:32-41.
 141. Weil ZM, Gaier KR, Karelina K (2014) Injury timing alters metabolic, inflammatory and functional outcomes following repeated mild traumatic brain injury. **Neurobiol. Dis.** 70:108-16.

**For more references, please visit our website
at www.fdneurotech.com.**

Notes

Notes
