



Quality & Excellence

FD Rapid GolgiStain™ Kit

A complete Golgi-Cox staining system for the study of the morphology of neurons and glia

User Manual
PK 401, Version 2009-01

FOR *IN VITRO* RESEARCH USE ONLY
not for diagnostic or other uses

FD NeuroTechnologies
Consulting & Services, Inc.

© 2002-2009 All Rights Reserved

Contents

I.	Introduction	4
II.	Kit Contents	5
III.	Materials Required but Not Included	5
IV.	Safety and Handling Precautions	6
V.	Tissue Preparation	6
VI.	Staining Procedure	11
VII.	References	13

I. Introduction

Golgi-Cox impregnation^{2,5} has been one of the most effective techniques for studying both the normal and abnormal morphology of neurons as well as glia. Using the Golgi technique, subtle morphological alterations in neuronal dendrites and dendritic spines have been discovered in the brains of animals treated with drugs as well as in the postmortem brains of patients with neurological diseases^{1,3}. However, the reliability and time-consuming process of Golgi staining have been major obstacles to the widespread application of this technique.

FD Rapid GolgiStain™ kit is designed based on the principle of the methods described by Ramón-Moliner⁵ and Glaser and Van der Loos⁴. This kit has not only dramatically improved and simplified the Golgi-Cox technique but has also proven to be extremely reliable and sensitive for demonstrating morphological details of neurons and glia, especially dendritic spines. The FD Rapid GolgiStain™ kit has been tested extensively on the brains from several species of animals as well as on the specimens of postmortem human brains (for photo samples and references using this kit, please visit our web site at www.fdtype.com).

II. Kit Contents

Store at room temperature

Solution A	250 ml
Solution B	250 ml
Solution C	250 ml x 2
Solution D	250 ml
Solution E	250 ml
Glass specimen retriever	2
Natural hair paintbrush	2
Dropping bottle	1
Plastic forceps	1
User manual	1

III. Materials Required but Not Included

1. Double distilled or Milli-Q water
2. Plastic/glass tubes or vials
3. Histological supplies and equipment:
 - Gelatin-coated microscope slides
 - Coverslips
 - Staining jars
 - Ethanol
 - Xylene
 - Permout[®]
 - A light microscope.

Permout[®] is a registered trademark of Fisher Scientific.

IV. Safety and Handling Precautions*

1. FD Rapid GolgiStain™ kit is made for *in vitro* research use only and not for drug, diagnostic or other uses.
2. The kit contains reagents that are toxic and harmful in contact with skin or by inhalation and may be fatal if ingested. Do not pipette by mouth. Avoid inhalation and contact with skin and eyes. In case of contact, wash immediately with generous amounts of water and seek medical advice. If swallowed, wash out mouth with water and immediately call a physician.
3. Perform experiment under a chemical hood. **Wear suitable protective clothing, gloves and eye/face protection while handling kit reagents.** Wash hands thoroughly after performing the experiment.

*Material safety data sheet (MSDS) is available at www.fdtypeurotech.com.

V. Tissue Preparation

The following instructions must be read before using this kit.

- All containers (plastic preferred) to be used should be cleansed and rinsed with distilled water.
- Do not use metal implements whenever Solutions A and B are present.
- Keep containers tightly closed at all times.
- Tissues treated with Solutions A and B, including sections should be protected from light whenever possible.

- The following procedure should be performed at room temperature unless specifically indicated.

 **Note**

FD Rapid GolgiStain™ kit has been proven to produce the best results in animal and postmortem human brain tissues prepared according to the following procedures. However, it may be used in tissue specimens prepared differently (please contact FD NeuroTechnologies for more information before use).

1. Experimental animals should be deeply anesthetized before killing. The animal brain (or postmortem human specimens) should be removed from the skull as quickly as possible but handled carefully to avoid damage or pressing of the tissue.

 **Note**

- *Perfusion of animals is not necessary. However, if you do perfuse animals (with paraformaldehyde only), tissue must not be postfixated. Important!*
 - *Large brain specimens should be sliced with a sharp blade into blocks of approximately 10 mm thickness.*
2. Rinse tissue briefly in double distilled or Milli-Q water to remove blood from the surface.
 3. Immerse tissue in the impregnation solution, made by mixing equal volumes of Solutions A and B, and store at room temperature for 2 weeks in the dark*. Replace the impregnation solution after the first 6 hours of immersion or on the next day.

* A 2-week impregnation time is satisfactory in most cases. However, variations in type and actual size of tissue may require a shorter or longer duration of impregnation to obtain the best results. The ideal time should be obtained by trial for each type of tissue, but 3 weeks is often sufficient for most tissues. Note that prolonging the impregnation time may increase background staining.

Note

- For the best results, the mixture of Solution A and B should be prepared at least 24 hours prior to use and left unstirred.
- It is important to use the top part of solution that is free of precipitate.
- The impregnation solution may be stored at room temperature for up to 1 month in the dark before use.
- Use at least 5 ml of the impregnation solution for each cubic cm of the tissue to be studied (i.e., the volume of the impregnation solution should be at least five times that of the tissue). Note that use a lesser volume of impregnation solution may decrease the sensitivity and reliability of staining.

Warning:

Solutions A and B (containing mercuric chloride, potassium dichromate and potassium chromate) are toxic in contact with skin and may be fatal if swallowed. The experiment should be performed under a chemical hood. Wear suitable protective clothing, gloves and eye/face protection while handling the reagents. **DO NOT POUR THE WASTE OF SOLUTIONS A AND B INTO THE SINK.** Collect the waste of these solutions in a bottle and

call your safety office or a licensed professional waste disposal service to dispose of this material.

4. Transfer tissue into Solution C and store at 4°C for at least 48 hours (up to 1 week) in the dark. Replace the solution after the first 24 hours of immersion or on the next day.
5. Sections (80 to 240 μm in thickness) may be cut slowly on a cryostat at -22°C (other types of microtomes, including sliding microtome and vibratome, can also be used, but contact FD Neuro-Technologies for more information before use). Each section should then be transferred with a glass specimen retriever (provided) and mounted on gelatin-coated microscope slides with Solution C (a dropping bottle is provided for easy dropping of Solution C onto the slides). Excess solution left on slide should be suctioned with a Pasteur pipette and then absorbed with a strip of filter paper (solution on slide must be wiped away as much as possible or sections will fall off slides). Allow sections to dry naturally at room temperature (do not use a fan or hot plate). Dried sections should be processed as soon as possible but may be stored in a slide box at room temperature for up to a week in the dark.

 **Note**

- *To prevent tissue from ice crystal damage and to preserve the best possible cell morphology, tissue should be frozen rapidly before sectioning with a cryostat. For*

example, tissue may be rapidly frozen as described below: place tissue in a plastic spoon and slowly dip into iso-pentane pre-cooled with dry ice (for the best results, temperature should be kept below -70°C and the dipping should take about 1 min, the slower the better). After the tissue is completely immersed in iso-pentane, keep it in iso-pentane for a few seconds and then place it on dry ice for another minute to ensure that the tissue is well frozen. Do not let tissue thaw before sections are cut.

- *The types of cryostat may vary, but all types should be able to cut thick sections (e.g. $100\ \mu\text{m}$). Contact FD NeuroTechnologies for technical assistance.*
- *If the cryostat has only one temperature control, set the cryostat temperature to -22°C at least 4 hours before cutting. If the cryostat has 2 temperature settings, set the chamber temperature 1 degree lower (colder) than that of the specimen head. Please note that -22°C is satisfactory in most cases. However, variations in type of cryostat and tissue may require a higher or lower chamber temperature in order to cut sections smoothly and without shattering.*
- *The following tips should be helpful for sectioning with a cryostat: 1) Mount tissue on specimen disc (also called chuck or stage) with distilled water, but make sure that the tissue is not thawed (may be done on dry ice). Tissue may also be mounted with any type of tissue freezing medium, including OCT, but avoid cutting through the medium. Do not embed the tissue in OCT. If the tissue has to be embedded for cutting, use TFM (TBS, Durham, NC, USA, Cat. #TFM-5) instead; 2) After the tissue is mounted on specimen disc, place the tissue on dry ice for 10 minutes, and then immediately install the speci-*

men disc together with the tissue on a cryostat and wait 5 minutes before cutting; 3) Cut a few sections slowly (do not use the anti-roll plate). If the tissue is too cold, e.g. sections show cracks which are parallel to the blade, wait a few more minutes before the next trial. Otherwise continue to cut.

VI. Staining Procedure

Do not let sections dry out between any steps during the staining and before coverslipping.

1. Rinse sections in double distilled or Milli-Q water 2 times, 2 minutes each.
2. Place sections in a mixture consisting of 1 part **Solution D**, 1 part **Solution E** and 2 parts double distilled or Milli-Q water for 10 minutes.

e.g. Solution D	10 ml
Solution E	10 ml
Double distilled water	20 ml

Note

- *The working solution should be prepared just before use and may be used for up to 200 sections (e.g. mouse brain) per 100 ml, depending on the size of sections.*
- *The bottle and staining jar containing the working solution must be covered to prevent vaporization of the reagent.*

- *For the best results, the working solution should be stirred frequently during incubation.*
3. Rinse sections in double distilled or Milli-Q water 2 times, 4 minutes each (distilled water should be renewed frequently).
 4. Counterstain sections with cresyl violet or thionin (optional step).

 **Note**

- *For the best results, step 3 must be prolonged, e.g. instead of 2 times and 4 minutes each, wash 4 times and 5 minutes each or longer.*
5. Dehydrate sections in 50%, 75% and 95% ethanol, 4 minutes each (do not skip any step).
 6. Dehydrate sections in absolute ethanol, 4 times, 4 minutes each (do not prolong).
 7. Clear in xylene, 3 times, 4 minutes each, and coverslip with PermOUNT[®].

 **Note**

- *Sections may be temporally stored in xylene for a few hours before coverslipping.*
- *For the best results, use undiluted PermOUNT[®].*
- *Golgi-stained sections should be protected from light whenever possible.*

VII. References

1. Robinson TE and Kolb B. Persistent structural modification in nucleus accumbens and prefrontal cortex neurons produced by previous experience with amphetamine. *J. Neurosci.* 17:8491-8497, 1997.
2. Corsi P. Camillo Golgi's morphological approach to neuroanatomy. In Masland RL, Portera-Sanchez A and Toffano G (eds.), *Neuroplasticity: a new therapeutic tool in the CNS pathology*, pp 1-7. Berlin: Springer, 1987.
3. Graveland GA, Williams RS and DiFiglia M. Evidence for degenerative and regenerative changes in neostriatal spiny neurons in Huntington's disease. *Science* 227:770-773, 1985.
4. Glaser ME and Van der Loos H. Analysis of thick brain sections by obverse-reverse computer microscopy: application of a new, high clarity Golgi-Nissl stain. *J. Neurosci. Methods* 4:117-125, 1981.
5. Ramón-Moliner E. The Golgi-Cox technique. In Nauta WJH and Ebesson SOE (eds.), *Contemporary Methods in Neuroanatomy*. pp 32-55, New York: Springer, 1970.

References using this kit:

1. Kurz JE, Moore BJ, Henderson SC, Campbell JN, Churn SB. A cellular mechanism for dendritic spine loss in the pilocarpine model of status epilepticus. ***Epilepsia*** 49(10):1696-710, 2008.
2. Chen Q, Xiong X, Lee TH, Liu Y, Wetsel WC, Zhang X. Neural plasticity and addiction: integrin-linked kinase and cocaine behavioral sensitization. ***J Neurochem.*** 2008 Sep 18. [Epub ahead of print]
3. Yamasaki N, Maekawa M, Kobayashi K, Kajii Y, Maeda J, Soma M, Takao K, Tanda K, Ohira K, Toyama K, Kanzaki K, Fukunaga K, Sudo Y, Ichinose H, Ikeda M, Iwata N, Ozaki N, Suzuki H, Higuchi M, Suhara T, Yuasa S, Miyakawa T. Alpha-CaMKII deficiency causes immature dentate gyrus, a novel candidate endophenotype of psychiatric disorders. ***Mol Brain*** 1(1):6, 2008.
4. Chen Q, Xiong X, Lee TH, Liu Y, Sun QA, Wetsel W, Zhang X. Adeno-associated virus-mediated ILK gene silencing in the rat NAc core. ***J Neurosci Methods*** 173(2):208-14, 2008.
5. Melendez-Ferro M, Perez-Costas E, Roberts RC. A new use for long-term frozen brain tissue: Golgi impregnation. ***J Neurosci Methods*** 2008 Aug 26. [Epub ahead of print]
6. Li CQ, Liu D, Huang L, Wang H, Zhang JY, Luo XG. Cytosine arabinoside treatment impairs the remote spatial memory function and induces dendritic retraction in the anterior cingulate cortex of rats. ***Brain Res Bull.*** 2008 Aug 25. [Epub ahead of print]
7. Zaja-Milatovic S, Gupta RC, Aschner M, Montine TJ, Milatovic D. Pharmacologic suppression of oxidative damage and dendritic degeneration following kainic acid-induced excitotoxicity in mouse cerebrum. ***Neurotoxicology*** 29(4):621-7, 2008.

8. An JJ, Gharami K, Liao GY, Woo NH, Lau AG, Vanevski F, Torre ER, Jones KR, Feng Y, Lu B, Xu B. Distinct role of long 3' UTR BDNF mRNA in spine morphology and synaptic plasticity in hippocampal neurons. **Cell** 134(1):175-87, 2008.
9. Samuels IS, Karlo JC, Faruzzi AN, Pickering K, Herrup K, Sweatt JD, Saitta SC, Landreth GE. Deletion of ERK2 mitogen-activated protein kinase identifies its key roles in cortical neurogenesis and cognitive function. **J Neurosci.** 28(27):6983-95, 2008.
10. Leão RN, Leão RM, da Costa LF, Rock Levinson S, Walmsley B. A novel role for MNTB neuron dendrites in regulating action potential amplitude and cell excitability during repetitive firing. **Eur J Neurosci.** 27(12):3095-108, 2008.
11. Liu QA, Shio H. Mitochondrial morphogenesis, dendrite development, and synapse formation in cerebellum require both Bcl-w and the glutamate receptor delta2. **PLoS Genet.** 4(6):e1000097, 2008.
12. Zhang L, Hernández VS, Medina-Pizarro M, Valle-Leija P, Vega-González A, Morales T. Maternal hyperthyroidism in rats impairs stress coping of adult offspring. **J Neurosci Res.** 86(6):1306-15, 2008.
13. Davis TH, Cuellar TL, Koch SM, Barker AJ, Harfe BD, McManus MT, Ullian EM. Conditional loss of Dicer disrupts cellular and tissue morphogenesis in the cortex and hippocampus. **J Neurosci.** 28(17):4322-30, 2008.
14. Savonenko AV, Melnikova T, Laird FM, Stewart KA, Price DL, Wong PC. Alteration of BACE1-dependent NRG1/ErbB4 signaling and schizophrenia-like phenotypes in BACE1-null mice. **Proc Natl Acad Sci U S A.** 105(14):5585-90, 2008.
15. Li X, Zhang P, Yang Y, Xiong Y, Qi Y, Hu H. Differentiation and developmental origin of cerebellar granule neuron ectopia in protein O-mannose UDP-N-acetylglucosaminyl transferase 1 knockout mice. **Neuroscience** 152(2):391-406, 2008.
16. Mishra A, Knerr B, Paixão S, Kramer ER, Klein R. The protein dendrite arborization and synapse maturation 1 (Dasm-1) is dispensable for dendrite arborization. **Mol Cell Biol.** 28(8):2782-91, 2008.
17. Quach TT, Massicotte G, Belin MF, Honnorat J, Glasper ER, Devries AC, Jakeman LB, Baudry M, Duchemin AM, Kolattukudy PE. CRMP3 is required for hippocampal CA1 dendritic organization and plasticity. **FASEB J.** 22(2):401-9, 2008.
18. Brun VH, Leutgeb S, Wu HQ, Schwarcz R, Witter MP, Moser EI, Moser MB. Impaired spatial representation in CA1 after lesion of direct input from entorhinal cortex. **Neuron** 57(2):290-302, 2008.
19. Fraser MM, Bayazitov IT, Zakharenko SS, Baker SJ. Phosphatase and tensin homolog, deleted on chromosome 10 deficiency in brain causes defects in synaptic structure, transmission and plasticity, and myelination abnormalities. **Neuroscience** 151(2):476-88, 2008.
20. Ye H, Tan YL, Ponniah S, Takeda Y, Wang SQ, Schachner M, Watanabe K, Pallen CJ, Xiao ZC. Neural recognition molecules CHL1 and NB-3 regulate apical dendrite orientation in the neocortex via PTP alpha. **EMBO J.** 27(1):188-200, 2008.
21. Dindot SV, Antalffy BA, Bhattacharjee MB, Beaudet AL. The Angelman syndrome ubiquitin ligase localizes to the synapse and nucleus, and maternal deficiency results in abnormal dendritic spine morphology. **Hum Mol**

- Genet.** 17(1):111-8, 2008.
22. Kikkawa YS, Pawlowski KS. Cochlear neuronal tracing for frequency mapping with DiI, NeuroVue, and Golgi methods. **Acta Otolaryngol Suppl.** 559:19-23, 2007.
 23. Kimura T, Yamashita S, Fukuda T, Park JM, Murayama M, Mizoroki T, Yoshiike Y, Sahara N, Takashima A. Hyperphosphorylated tau in parahippocampal cortex impairs place learning in aged mice expressing wild-type human tau. **EMBO J.** 26(24):5143-52, 2007.
 24. Gunnarsen JM, Kim MH, Fuller SJ, De Silva M, Britto JM, Hammond VE, Davies PJ, Petrou S, Faber ES, Sah P, Tan SS. Seiz-6 proteins affect dendritic arborization patterns and excitability of cortical pyramidal neurons. **Neuron** 56(4):621-39, 2007.
 25. Wallace M, Frankfurt M, Arellano A, Inagaki T, Luine V. Impaired recognition memory and decreased prefrontal cortex spine density in aged female rats. **Ann N Y Acad Sci.** 1097:54-7, 2007.
 26. Ampuero E, Dagnino-Subiabre A, Sandoval R, Zepeda-carreno R, Sandoval S, Viedma A, Aboitiz F, Orrego F and Wyneken U. Status epilepticus induces region-specific changes in dendritic spines, dendritic length and TrkB protein content of rat brain cortex. **Brain Research** 1150:225-238, 2007.
 27. Conrad CD, McLaughlin KJ, Harman JS, Foltz C, Wieczorek L, Lightner E and Wright RL. Chronic glucocorticoids increase hippocampal vulnerability to neurotoxicity under conditions that produce CA3 dendritic retraction but fail to impair spatial recognition memory. **J. Neuroscience** 27:8278-8285, 2007.
 28. Cornejo BJ, Mesches MH, Coultrap S, Browning MD and Benke TA. A single episode of neonatal seizures permanently alters glutamatergic synapses. **Annals of Neurology** 61:411-426, 2007.
 29. Fünfschilling U, Saher G, Xiao L, Möbius W and Nave A. Survival of adult neurons lacking cholesterol synthesis in vivo. **BMC Neuroscience** 8:1, 2007.
 30. Johnson DA, Zhang J, Frase S, Wilson M, Rodriguez-Galindo C and Dyer MA. Neuronal differentiation and synaptogenesis in retinoblastoma. **Cancer Research** 67:2701-2711, 2007.
 31. Keys DA, Tian G, Poirier K, Huang GJ, Siebold C, Cleak J, Oliver PL, Fray M, Harvey RJ, Molnar Z, Pinon MC, Dear N, Valdar W, Brown SDM, Davies KE, Rawlins JNP, Cowan NJ, Nolan P, Chelly J and Flint J. Mutations in α -tubulin cause abnormal neuronal migration in mice and lissencephaly in humans. **Cell** 128:45-57, 2007.
 32. Kennedy NJ, Martin G, Ehrhardt AG, Cavanagh-Kyros J, Kuan CY, Rakic P, Flavell RA, Treisman SN and Davis RJ. Requirement of JIP scaffold proteins for NMDA-mediated signal transduction. **Genes & Development** 21:2336-2346, 2007.
 33. Li L, Yun SH, Keblesh J, Trommer BL, Xiong H, Radulovic J and Tourtellotte WG. Egr3, a synaptic activity regulated transcription factor that is essential for learning and memory. **Molecular and Cellular Neuroscience** 35:76-88, 2007.
 34. McLaughlin KJ, Gomez JL, Baran SE and Conrad CD. The effects of chronic stress on hippocampal morphology and function: an evaluation of chronic restraint paradigms. **Brain Res.** 1161:56-64, 2007.
 35. Nielsen JV, Nielsen FH, Ismail R, Noraberg J and Jensen NA. Hippocampus

- like corticoneurogenesis induced by two isoforms of the BTB-zinc finger gene *Zbtb20* in mice. **Development** 134:1133-1140, 2007.
36. Ostwald D, Karpac J, and Hochgeschwender U. Effects on hippocampus of lifelong absence of glucocorticoids in the pro-opiomelanocortin null mutant mouse reveal complex relationship between glucocorticoids and hippocampal structure and function. **Journal of Molecular Neuroscience** 28:291-302, 2007.
 37. Quach TT, Massicotte G, Belin MF, Honnorat J, Glasper ER, Devries AC, Jakeman LB, Baudry M, Duchemin AM and Kolattukudy PE. CRMP3 is required for hippocampal CA1 dendritic organization and plasticity. **The FASEB J.** doi: 10.1096/fj.07-9012com, 2007.
 38. Sadakata T, Kakegawa W, Mizoguchi A, Washida M, Katoh-Semba R, Shutoh F, Okamoto T, Nakashima H, Kimura K, Tanaka M, Sekine Y, Ito-hara S, Yuzaki M, Nagao S and Furuich T. Impaired cerebellar development and function in mice lacking CAPS2, a protein involved in neurotrophin release. **J. Neuroscience** 27:2472-2482, 2007.
 39. Sarti F, Borgland SL, Kharazia VN and Bonci A. Acute cocaine exposure alters spine density and long-term potentiation in the ventral tegmental area. **Eur J. Neuroscience** 26:749-756, 2007.
 40. Shu Y, Duque A, Yu Y, Haider B and McCormick DA. Properties of action potential initiation in neocortical pyramidal cells: evidence from whole cell axon recordings. **J. Neurophysiology** 97:746-760, 2007.
 41. Soderling SH, Guire ES, Kaech S, White J, Zhang F, Schutz K, Langeberg LK, Banker G, Raber J and Scott JD. A WAVE-1 and WRP signaling complex regulates spine density, synaptic plasticity, and memory. **J. Neuroscience** 27:355-365, 2007.
 42. Sun M, Thomas MJ, Herder R, Bofenkamp ML, Selleck SB and O'Connor MB. Presynaptic contributions of chordin to hippocampal plasticity and spatial learning. **J. Neuroscience** 27:7740-7750, 2007.
 43. Wallace M, Frankfurt M, Arellanos A, Inagaki T and Luine V. Impaired recognition memory and decreased prefrontal cortex spine density in aged female rats. **Ann. N. Y. Acad. Sci.** 1097: 54-57, 2007.
 44. Watase K, Gatchel JR, Sun Y, Emamian E, Atkinson R, Richman R, Mizusawa H, Orr HT, Shaw C and Zoghbi HY. Lithium therapy improves neurological function and hippocampal dendritic arborization in a spinocerebellar ataxia type 1 mouse model. **PLoS Medicine** 4:836-847, 2007.
 45. Adhami F, Liao GH, Morozov YM, Schloemer A, Schmthorst VJ, Lorenz JN, Dunn RS, Vorhees CV, Wills-Karp M, Degen JL, Davis RJ, Mizushima N, Rakic P, Dardzinski BJ, Holland SK, Sharp FR and Kuan CY. Cerebral ischemia-hypoxia induces intravascular coagulation and autophagy. **American J. Pathology** 169:566-583, 2006.
 46. Buss RR, Gould TW, Ma J, Vinsant S, Prevette D, Winseck A, Toops KA, Hammarback JA, Smith TL and Oppenheim RW. Neuromuscular development in the absence of programmed cell death: phenotypic alteration of motoneurons and muscle. **J. Neuroscience** 26:13413-13427, 2006.
 47. Dang MT, Yokoi FY, Yin HH, Lovinger DM, Wang Y and Li Y. Disrupted motor learning and longterm synaptic plasticity in mice lacking NMDAR1 in the striatum. **Proc. Natl. Acad. Sci. USA** 103:15254-15259, 2006.
 48. Elia LP, Yamamoto M, Zang K and Reichardt LF. p120 catenin regulates dendritic spine and synapse development through Rho-family GTPases and

- cadherins. **Neuron** 51:43-56, 2006.
49. Huang Z, Shimazu K, Woo NH, Zang K, Müller U, Lu B and Reichardt LF. Distinct roles of the β 1-class integrins at the developing and the mature hippocampal excitatory synapse. **J. Neuroscience** 26:11208-11219, 2006.
 50. Inan M, Lu HC, Albright MJ, She WC and Crair MC. Barrel map development relies on protein kinase A regulatory subunit II β -mediated cAMP signaling. **J. Neuroscience** 26:4338-4349, 2006.
 51. Jacobsen JS, Wu CC, Redwine JM, Comery TA, Arias R, Bowlby M, Martone R, Morrison JH, Pangalos, MN, Reinhart PH and Bloom FE. Early-onset behavioral and synaptic deficits in a mouse model of Alzheimer's disease. **Proc. Natl. Acad. Sci. USA** 103:5161-5166, 2006.
 52. Kleen JK, Sitomer MT, Killeen PR and Conrad CD. Chronic stress impairs spatial memory and motivation for reward without disrupting motor ability and motivation to explore. **Behavioral Neuroscience** 120:842-851, 2006.
 53. Luine V, Attalla S, Mohan G, Costa A and Frankfurt M. Dietary phytoestrogens enhance spatial memory and spine density in the hippocampus and prefrontal cortex of ovariectomized rats. **Brain Research** 1126:183-187, 2006.
 54. Maezawa I, Zaja-Milatovic S, Milatovic D, Stephen C, Sokal I, Maeda N, Montine TJ and Montine KS. Apolipoprotein E isoform-dependent dendritic recovery of hippocampal neurons following activation of innate immunity. **J. Neuroinflammation** 3:21, 2006.
 55. Melendez-Ferro M, Perez-Costas E and Roberts RC. Golgi staining of long-term stored human brain tissue. **FENS Abstr.** 3:A059.20, 2006.
 56. Moretti P, Levenson JM, Battaglia F, Atkinson R, Teague R, Antalffy B, Armstrong D, Arancia O, Sweatt JD and Zoghbi HY. Learning and memory and synaptic plasticity are impaired in a mouse model of Rett syndrome. **J. Neuroscience** 26:319-327, 2006.
 57. Poet M, Kornak U, Schweizer M, Zdebek AA, Scheel O, Hoelter S, Wurst W, Schmitt A, Fuhrmann JC, Planells-Cases R, Mole SE, Hubner CA and Jentsch TJ. Lysosomal storage disease upon disruption of the neuronal chloride transport protein CIC-6. **Proc. Natl. Acad. Sci. USA** 103:13854-13859, 2006.
 58. Wallace M, Luine V, Arellanos A and Frankfurt M. Ovariectomized rats show decreased recognition memory and spine density in the hippocampus and prefrontal cortex. **Brain Research** 1126:176-182, 2006.
 59. Watson KK, Jones TK and Allman JM. Dendritic architecture of the von Economo neurons. **Neuroscience** 141:1107-1112, 2006.
 60. Björkblom B, Östman N, Hongisto V, Komarovski V, Filén JJ, Nyman TA, Kallunki T, Courtney MJ and Coffey ET. Constitutively active cytoplasmic c-jun N-terminal kinase 1 is a dominant regulator of dendritic architecture: role of microtubule-associated protein 2 as an effector. **J. Neuroscience** 25:6350-6361, 2005.
 61. Dagnino-Subiabre A, Terreros G, Carmona-Fontaine C, Zepeda R, Orellana JA, Diaz-Veliz G, Mora S and Aboitiz F. Chronic stress impairs acoustic conditioning more than visual conditioning in rats: morphological and behavioural evidence. **Neuroscience** 135:1067-1074, 2005.
 62. Grillet N, Pattyn A, Contet C, Kieffer BL, Goridis C and JF Brunet. Generation and characterization of Rgs4 mutant mice. **Molecular and Cellular**

- Biology 25:4221-4228, 2005.
63. Gu X, Li C, Wei W, Lo V, Gong S, Li SH, Iwasato T, Itohara S, Li XJ, Mody I, Heintz N and Yang XW. Pathological cell-cell interactions elicited by a neuropathogenic form of mutant huntingtin contribute to cortical pathogenesis in HD mice. **Neuron** 46:433-444, 2005.
 64. Ishikura N, Clever JL, Bouzamondo-Bernstein E, Samayoa E, Prusiner SB, Huang EJ and DeArmond SJ. Notch-1 activation and dendritic atrophy in prion disease. **Proc. Natl. Acad. Sci. USA** 102:886-891, 2005.
 65. Laub F, Lei L, Sumiyoshi H, Kajimura D, Dragomir C, Smaldone S, Puche AC, Petros TJ, Mason C, Parada LF and Ramirez F. Transcription factor KLF7 is important for neuronal morphogenesis in selected regions of the nervous system. **Molecular and Cellular Biology** 25:5699-5711, 2005.
 66. McLaughlin KJ, Baran SE, Wright RL and Conrad CD. Chronic stress enhances spatial memory in ovariectomized female rats despite CA3 dendritic retraction: possible involvement of CA1 neurons. **Neuroscience** 135:1045-1054, 2005.
 67. Niewmierzycka A, Mills J, St-Arnaud R, Dedhar S and Reichardt LF. Integrin-linked kinase deletion from mouse cortex results in cortical lamination defects resembling cobblestone lissencephaly. **J. Neuroscience** 25:7022-7031, 2005.
 68. Pyter LM, Reader BF and Nelson RJ. Short photoperiods impair spatial learning and alter hippocampal dendritic morphology in adult male white-footed mice (*Peromyscus leucopus*). **J. Neuroscience** 25:4521-4526, 2005.
 69. Ramanan N, Shen Y, Sarsfield S, Lemberger T, Schütz G, Linden DJ and Ginty DD. SRF mediates activity-induced gene expression and synaptic plasticity but not neuronal viability. **Nature Neuroscience** 8:759-767, 2005.
 70. Ren-Patterson RF, Cochran LW, Holmes A, Sherrill S, Huang SJ, Tolliver T, Lesch KP, Lu B and Murphy DL. Loss of brain-derived neurotrophic factor gene allele exacerbates brain monoamine deficiencies and increases stress abnormalities of serotonin transporter knockout mice. **J. Neurosci. Res.** 79:756-771, 2005.
 71. Amateau SK and McCarthy MM. Induction of PGE₂ by estradiol mediates developmental masculinization of sex behavior. **Nature Neuroscience** 7:643-650, 2004.
 72. Beggs HE, Schahin-Reed D, Zang K, Goebels S, Nave KA, Gorski J, Jones KR, Sretavan D and Reichardt LF. FAK deficiency in cells contributing to the basal lamina results in cortical abnormalities resembling congenital muscular dystrophies. **Neuron** 40:501-514, 2003.
 73. Dahl JP, Wang-Dunlop J, Gonzales C, Goad MEP, Mark RJ and Kwak SP. Characterization of the WAVE1 knock-out mouse: implications for CNS development. **J. Neuroscience** 23:3343-3352, 2003.
 74. Mlatovic D, Zaja-Milatovic S, Montine KS, Horner PJ, and Montine TJ. Pharmacologic suppression of neuronal oxidative damage and dendritic degeneration following direct activation of glial innate immunity in mouse cerebrum. **J. Neurochem.** 87:1518-1526, 2003.

**For more references, please visit our website at
www.fднеurotech.com.**

Notes

Notes
